

GARIS PANDUAN PENGHANTARAN SPESIMEN KE SEKSYEN KAWALAN MUTU

Makmal	Ujian	Jenis Spesimen	Bekas spesimen	Isipadu (ml)	Borang	TAT	Syarat Khusus
Hematologi	Full Blood Count (FBC) – pre and post apheresis donor	Whole blood	K ₂ EDTA	2ml – 4ml	QC for Blood Products (PDN/QA/QP-01/01)	Urgent 30 min Normal request 3 hour	Tiada
	QC for blood and blood component						
	Full Blood Count (FBC)	Whole blood, Red Cell, Platelet	K ₂ EDTA	2 ml – 4 ml	QC for Blood Products (PDN/QA/QP-01/01)	24 jam	Tiada
	Factor VIII	Plasma	Plain tube	2 ml – 4 ml		24 jam	Tiada
	Fibrinogen	Plasma	Plain tube	2 ml – 4 ml		24 jam	Tiada
	Leukocyte Count	Whole blood, Red Cell	Plain tube	2 ml – 4 ml		24 jam	Tiada
	pH	Platelet	Plain tube	2 ml		30 minit	Tiada
	% Hemolysis	Whole Blood, Red Cell	Plain tube	4 ml		1 jam	Tiada
	Osmolality	Deglycerolized product	Plain tube	4 ml		1 jam	Tiada
	Sterility	All products	Blood bag/ transfer bag	20 ml		7 hari	Steril